

TALLADIUM U.K.



37B New Road, Langley, Berks SL3 8JJ
Tel: 01753 582 788 Fax: 01753 582 422
info@talladiumuk.co.uk www.talladiumuk.co.uk

idea
Inspiring Dental
Education & Advice

idea

Learning Provider Partner

FLEXIBLE DENTURE COURSE

6 hours verifiable CPD

A ONE DAY COURSE THAT WILL COVER:

- *The indications and contra-indications for flexible partials*
- *A guide to the different materials & the various appliances that can be made*
- *The outline of the procedure for all these materials from design to finish including:*
 - *Blocking undercuts*
 - *Model duplication*
 - *Teeth set up*
 - *Wax up*
 - *Investing and boiling out*
 - *Injecting*
 - *Finishing and polishing (including a guide to the best tools for the job)*
- *Common mistakes and pitfalls*
- *Practical tips and techniques at each stage from your tutor, Mr Slavko Kos RDT MDTA, a dental technician with over 40 years experience, who has been using flexible materials for over 10 years.*

ON COMPLETION OF THE COURSE YOU WILL RECEIVE:

- *A demonstration model to show your dentists*
- *Promotional literature for dentists for all materials demonstrated on the course*
- *Manuals and instructions for all materials demonstrated on the course*
- *IDEA-issued certificate*
- *6 hours verified CPD*

COURSE DETAILS:

Location: Langley, Berkshire
Cost: £350
Course dates: Several times a month including Saturdays

CALL NOW ON 01753 582 788 TO RESERVE YOUR PLACE!



Don't be tied to just one flexible material - the **J-100 Evolution** can process them all !



Combination with tooth colour clasps



Flexident

From the manufacturers of the Talladium Tilite alloy range, comes

FLEXIDENT - for lightweight, strong, durable flexible partials and dentures.


Available as single cartridges.

Shades: light pink, regular pink, dark pink



Pressing Dental, the manufacturers of the J-100 Evolution Injection System, have created two new thermoplastic materials. Both are available in 250g packs.

FLEXI-J is a nylon-based thermoplastic. Type 1 is more elastic and translucent, whilst Type 2 is stiffer and totally



transparent. The two materials can be mixed to create an intermediary elasticity.

F.J.P is the result of years of research and experimentation with nylon-type prostheses. Created as an alternative compound to nylon, with improved features and advantages,

F.J.P is:
Moderately elastic
Semi-transparent
Precise
Colour stable
Easy to process
Easy to polish
Adhesion to acrylic possible



FLEXITE PLUS super flexible nylon for partials
Shades: light pink, translucent pink, dark pink, ethnic-purple, coe-ethnic moderate, ethnic-brown, natural (clear)

FLEXITE M.P. for full dentures & TMJs
Shades: pink, lucitone pink, coe brown ethnic, crystal clear and tooth colour shades

FLEXITE NORTHERM monomer and acrylic free for patients with allergies
Shades: pink, clear & ethnic



FLEXITE SUPREME can be used in place of chrome
Shades: natural (clear), light pink, dark pink, ethnic brown, coe 22 ethnic



FOR FURTHER INFORMATION OR TO BOOK A COURSE CALL **01753 582 788**

